OxyFerm™ DO Sensor

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Table of Contents

1	Introduction	4
2	Liability	4
3	Intended Use	4
4	Safety Precautions	5
5	Initial Operation	5
6	Electrical connection	6
7	Getting ready to take measurements	6
8	Dismantling the sensor	6
9	Sterilization	7
10	Sensor check and maintenance	7
11	Replacing the electrolyte and the membrane cartridge	8
12	Replacing the Cathode	9
13	Disposal	10
14	Instructions for use in potentially explosive atmospheres	10
	14.1 Assembly	10
	14.2 Connection to Transmitter	11
15	Accessories	12
16	Technical Data	12

Important Notice

Reproduction of any part of this manual in any form whatsoever without the express written consent of Hamilton Bonaduz AG is forbidden. The contents of this manual are subject to change without notice. Technical changes reserved.

All efforts have been made to ensure the accuracy of the contents of this manual. However, should any errors be detected, Hamilton Bonaduz AG would greatly appreciate being informed of them. The above notwithstanding, Hamilton Bonaduz AG can assume no responsibility for any errors in this manual or their consequences.

Copyright © 2016 Hamilton Bonaduz AG, Bonaduz Switzerland.
All rights reserved.

1 Introduction

These operating instructions refer to the sterilizable Hamilton oxygen sensors OxyFerm FDA, OxyFerm FDA VP and OxyFerm XL.

Name	Ref	ATEX
OxyFerm FDA 120	237 450	From serial number 4000
OxyFerm FDA 160	237 455	From serial number 4000
OxyFerm FDA 225	237 452	From serial number 4000
OxyFerm FDA 325	237 453	From serial number 4000
OxyFerm FDA 425	237 454	From serial number 4000
OxyFerm FDA VP 120	237 540	From serial number 4000
OxyFerm FDA VP 160	237 541	From serial number 4000
OxyFerm FDA VP 225	237 542	From serial number 4000
OxyFerm FDA VP 325	237 543	From serial number 4000
OxyFerm FDA VP 425	237 544	From serial number 4000
OxyFerm XL	237 175 - OP	From serial number 4000
OxyFerm XL 150	237 170	From serial number 4000
OxyFerm XL 300	237 174	From serial number 4000

Hamilton oxygen sensors are quality product produced in accordance with the latest scientific and technical knowledge. Follow the instructions given here and you can be assured of maximum safety and durability.

These instructions should be read, understood and followed by all staff using the device. Hamilton can assume no responsibility for damage and operational disruptions arising from failure to observe these instructions.

2 Liability

The liability of Hamilton Bonaduz AG is detailed in the document "General Terms and Conditions of Sale and Delivery (GTS)", chapter 12.

Hamilton is expressly not liable for direct or indirect losses arising from use of the sensors. It must in particular be ensured in this context that malfunctions can occur on account of the inherently limited useful life of sensors contingent upon their relevant applications. The user is responsible for the calibration, maintenance and regular replacement of the sensors. In the case of critical sensor applications, Hamilton recommends using back-up measuring points in order to avoid consequential damage. The user is responsible for taking suitable precautions in the event of a sensor failure.

3 Intended Use

OxyFerm sensors are designed to measure the partial pressure of dissolved oxygen in liquids. The main applications for this sensor are in biotechnology and pharmaceutical industries. The sensor's outstanding feature is its long-term stability even when sterilized or autoclaved often.

The OxyFerm sensor has a built-in temperature detector (NTC 22 kOhm). This temperature detector is to be used only for compensation of the oxygen signal, not for controlling the process temperature.

If the sensor is used in potentially explosive atmospheres, refer to the section "Instructions for use in potentially explosive atmospheres".

4 Safety Precautions

This sensor is to be used only as intended and is to be kept in a condition that ensures complete safety. The specifications given in the section "Technical Data" as regards temperature, pressure etc. may under no circumstances be exceeded. Inappropriate use or misuse can be dangerous.

Assembly and maintenance may only be done by trained personnel. Since the inner part of the sensor is glass, handle it with care. The top of the sensor can be damaged by impacts.

When changing the electrolyte, it is recommended to wear safety goggles and protective gloves. Contact of the electrolyte with skin or eyes is to be avoided (see Safety Data Sheet Ref 608914).

If mechanical damage to the sensor occurs, electrolyte may leak out.

Ensure that, when screwed into the process, the PG13,5 thread and the O-ring incur no damage. O-rings are consumable parts and must be replaced regularly – at least once a year. Even where all necessary safety measures have been complied with, there is still a further risk involving leaks or mechanical damage to the armature. Where there are seals or screws, gases or liquids could leak out undetected. Before dismantling the sensor, make sure that there is no over-pressure in the process.


$T_{\max} = 130^{\circ}\text{C}$

$P_{\max} = 4 \text{ bar}$

5 Initial Operation

During unpacking please check for possible mechanical defects. In case of a complaint please return the sensor in its original packing to your Hamilton partner.

1 Connector

3 Sensor Shaft

5 G1¼" Thread

2 PG13,5 Thread

4 Membrane cartridge

6 O-ring


⚠ ATTENTION! OxyFerm XL sensors have a 25 mm diameter process connection. Check that the O-ring (Pos. 6) is in the right position.

6 Electrical Connection

Please use the appropriate Hamilton connection cables (see Accessories).


OxyFerm FDA, OxyFerm XL

- A Cathode
- B Anode
- C, D NTC 22 kOhm

OxyFerm FDA VP

- A Cathode
- B Anode
- C Not connected
- D Shaft*
- E, F NTC 22 kOhm

* We recommend to connect Pin D to the system potential earth (PE)

7 Getting ready to take measurements

1. Remove the protective watering cap.
2. Unscrew the membrane cartridge and pre-rinse it with electrolyte. The OxyFerm sensors have no electrolyte in the membrane cartridge when they are new (exception OxyFerm FDA). Then fill it with 1.5 ml of electrolyte. Carefully screw the membrane cartridge tight again.
3. Mount the sensor by using the thread PG13,5 / G1¼. Do not use an upside down mounting. Hamilton has developed a special version for upside down mounting. Please contact Hamilton if a corresponding application is desired.
4. Connect the OxyFerm to the amplifier. Its polarization voltage should be -670 +/-50 mV. After a period of 2 hours minimum the sensor is ready for calibration. Such a polarization time is necessary to get stable signals. The following rule is valid if the sensor is disconnected from the amplifier for a short time: Polarization time = 2 x time of disconnecting, but no more than 2 hours.
5. The calibration is preferably performed in water vapour-saturated ambient air. When calibrated in normal ambient air only, a negligible measuring error occurs.
6. If the value is too small or the response sluggish, see "Sensor Check and Maintenance", point 2.

8 Dismantling the sensor

Unscrew the PG13,5 / G1¼ thread and pull out the sensor. Never turn the whole sensor. Make sure that there is no process pressure!


9 Sterilization

If the sensor is autoclaved, it is advantageous to protect the connector with a protective cap (Ref 237 137). You do not need a protective cap if using OxyFerm VP. No protection is necessary when performing an in-situ sterilization. If the connector is partially filled with water despite the protective cap, please rinse with ethanol and dry with a hair dryer. The drying is necessary to avoid corrosion or insulation problems.

The contacts have to be dry and clean before connecting to the cable.

10 Sensor check and maintenance

1. Calibrate the polarized sensor in air and adjust to 100%. Put the sensor in nitrogen gas (i.e. by using a bag filled with nitrogen). After 1 minute in nitrogen the sensor must show a value below 2% of the value in air. A too high current in nitrogen is due to an exhausted electrolyte or a deformed or defective gas-permeable membrane. In such a case replace the membrane cartridge as well as the electrolyte solution.
2. If there is a sluggish response and/or a too small current in air:
 - a) Go to "Replacing the electrolyte..." to replace the electrolyte only.
 - b) If still not calibrating well, go to "Replacing the electrolyte..." to replace the membrane cartridge as well as the electrolyte.
 - c) If the problem persists, this is an indication of a contaminated cathode. In this case the cathode and the front-end of the glass body need to be cleaned with the polishing strip included in the membrane kit. First, take the dry polishing strip in one hand; then press the sensor lightly into the polishing strip and turn it 10-20 times. Afterwards, rinse the glass body under running water and dry it carefully.


- d) If the problem still persists, you need to replace the cathode, see “Replacing the cathode”.

⚠ ATTENTION! Mechanical polishing of the cathode should only be done with Hamilton polishing cloths. There is a polishing cloth included in the membrane kit. Do not touch the anode wire. Do not bend the glass body.

11 Replacing the electrolyte and the membrane cartridge

The OxyFerm membrane cartridge is replaced as follows:

1. Put the sensor in an upright position and carefully unscrew the membrane cartridge.
2. Carefully clean the tip of the glass body with a soft tissue.


3. Check the small O-ring above the glass body. Replace if damaged.
4. Use the 1 ml plastic pipette of the membrane kit to pipette 1.5 ml of electrolyte (Oxylyte, Ref 237 118) into the new membrane cartridge.
5. Carefully screw the cartridge onto the sensor shaft. Any spillage of electrolyte should be rinsed away with water (see “Safety aspects”).

12 Replacing the Cathode

If the cathode is going to be changed, the membrane must also be changed to get the full performance of the sensor.

1. Put the sensor in an upright position and carefully unscrew the membrane cartridge.
2. Flush the replaceable cathode with deionized water, then dry the metal parts. Do not touch anode and cathode.
3. Hold the cathode on the metal part (B) in front of the thread (A) and the sensor on the shaft (C) and pull apart. Just pull, do not turn!
4. Please check that all contacts are dry and clean.
5. Replace cathode with a new one (Ref 237306 (blue) for OxyFerm; Ref 237427 (yellow) for OxyGold B; Ref 237437 (black) for OxyGold G). Rotate the cathode until the correct position is found and then push the cathode into the shaft.
6. Check the small O-ring (D) above the cathode. Replace if damaged.
7. Use the 1 ml plastic pipette of the membrane kit to pipette 1.5 ml of electrolyte into the new membrane cartridge. The pipette must NOT touch the membrane itself!
8. Screw the cartridge onto the sensor shaft. Any spillage of electrolyte should be rinsed away with water (see “Safety aspects”).

13 Disposal

The design of Hamilton sensors minimizes environmental impact. According to the EU directive 2002/96/EC the Hamilton sensors should be disposed of as waste of electrical and electronic equipment, and not in municipal waste, or it can be sent back to Hamilton for disposal.

14 Instructions for use in potentially explosive atmospheres

ATEX / IECEx marking:

Gas: CE 0035 ⚠ II 1/2 G Ex ia IIC T4/T5/T6 Ga/Gb

Dust: CE 0035 ⚠ II 1/2 D Ex ia IIIC T x °C Da/Db

Manufacturer: Hamilton Bonaduz AG, 7402 Bonaduz, Switzerland

EC type examination report: TÜV 03 ATEX 7005 X

IECEx certificate of conformity: IECEx TUR 14.0001 X

EC type examination report and IECEx certificate of conformity can be downloaded from www.hamiltoncompany.com.

The conditions described in the ATEX / IECEx certificate must be respected.

⚠ ATTENTION! In case a gas atmosphere and a dust atmosphere are or could be present at the same time, the risk of explosion must be examined carefully and special precautions may be necessary.

14.1 Assembly

- 1) The operator of equipment in potentially explosive atmospheres is responsible for ensuring that all components of the system are certified for that area classification and are compatible with each other.
- 2) In gas atmospheres, the Ex approval is not dependent on any conditions regarding mounting of the sensor. In dust atmospheres, however, special restrictions may apply when using armatures made of plastic material. Use with particular Hamilton armatures can lead to an increase in the allowable power of the transmitter. See section „Connection to Transmitter“ below.
- 3) O-rings having a sealing function between Ex-zone 0 and 1 or Ex-zone 20 and 21, respectively and must be replaced at each dismantling of the sensor.
- 4) Sensors, transmitter and equipment are required to be set up within one balanced potential system.
- 5) When the sensor is assembled, the ATEX / IECEx sticker is to be attached to the cable in an easily visible place, as near as possible to the sensor itself. This sticker is to indicate that an ATEX / IECEx approved sensor is in use. The sticker should not be removed at any time.

14.2 Connection to Transmitter

OxyFerm electrodes are suitable for connection to an intrinsically safe electrical circuit with protection level ia. The operator of the equipment must ensure that the allowable electrical values for the sensor all exceed those for the transmitter. None of the values cited for electrical power, voltage and current may be exceeded in sum (measurement and temperature circuit together). Electric values for OxyFerm are as follows:

U = 24 V; I = 173 mA; P: Function of process temperature. The maximum allowable electric power P of the transmitter (measurement and temperature circuit together) depends on the desired process temperature, and the coupling situation.

If OxyFerm is used with a RetractoFit or FlexiFit type armature, as per the table below, the values of type 5 are the correct ones. In all other coupling situations, especially for OxyFerm XL, the values of type 4 are the correct ones.

Coupling Situation for Type 5	OxyFerm 120	OxyFerm 225	OxyFerm XL
RetractoFit Series*	–	yes	–
RetractoFit Bio Series*	–	yes	–
FlexiFit Series*	yes	–	–

For gas atmospheres:

Type 4	P ≤ 30 mW	P ≤ 60 mW	P ≤ 90 mW
T4	129°C	129°C	128°C
T5	94°C	94°C	93°C
T6	68°C	57°C	46°C

Type 5	P ≤ 250 mW
T4	125°C
T5	90°C
T6	75°C

For dust atmospheres:

Table for the calculation of the maximum surface temperature „x“ of the sensor as a function of the maximum electrical power of the transmitter „P_i“ and ambient/process temperatures “T_a”. „x“ has to be smaller than the ignition temperatures of dust involved.

Type	P _i ≤ 30 mW	P _i ≤ 60 mW	P _i ≤ 90 mW
4	x = T _a + 12°C	x = T _a + 23°C	x = T _a + 34°C

Type	P _i ≤ 250 mW
5	x = T _a + 5°C

The temperature limit stated on the sensor (130°C) must always be respected.

* several reference numbers possible

15 Accessories

Name	Ref
Cable (1 m)	355 087
Cable VP 6 (1 m)	355 108
Cable (3 m)	355 088
Cable, VP 6 (3 m)	355 109
Cable (5 m)	355 089
Cable, VP 6 (5 m)	355 110
Cable, VP 6 (10 m)	355 111
Cable, VP 6 (20 m)	355 112
OxyFerm Membrane Kit	237 123
CIP Membrane Kit	237 126
FDA Membrane Kit	237 140
Oxylyte Electrolyte 50 ml	237 118
Autoclavation Cap	242 000
Replacement Cathode OxyFerm (blue)	237 306

16 Technical Data

Electrode system	silver-platinum combination
Electrolyte	Oxylyte
Membrane	Optiflow (FDA)
Temperature sensor	NTC 22 kOhm
Polarization voltage and time	-670 +/- 50 mV, ≥ 2 hours
Current in air (at 25 °C)	40 to 80 nA
Operating temperature range	0 to 130 °C
Storage temperature	-10 to 60 °C
Pressure range	0 to 4 bar
Diameter	12 mm standard; 25 mm XL-Type
a-length	several versions available
Process connection	PG13,5 - thread standard; G 1¼" XL-Type
Wetted parts	see certificate
Regenerating	membrane cartridge, electrolyte replaceable and cathode
Response time $t_{98\%}$	30 to 60 s at 25 °C, from air to nitrogen
Temperature response	3.1 %/K
Required Flow	≥ 0.03 m/s
Flow response	< 5% at 25 °C
Oxygen consumption	ca. 20 ng/h in air at 25 °C
Residual current	< 0.1% of current on air
Zero shift	< 0.1% of current on air
Drift at room temperature under constant conditions	< 1% per week

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Copyright © 2016 Hamilton Bonaduz AG, Bonaduz Schweiz.
© 2016 Hamilton Bonaduz AG. All rights reserved.